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Table 3-1: Typology of Liquid Biofuel

Fuel	Feedstock	Regions where currently produced	Greenhouse Gas (GHG) reduction impacts vs. petroleum fuel	Costs	Biofuel yield per hectare of land	Land Types
Ethanol	Grains (wheat, maize)	US, Europe, China	Low-moderate	Moderate	Moderate	Croplands
Ethanol	Sugar cane	Brazil, India, Thailand	High	Low	High	Croplands
Ethanol	Biomass (cellulose)	None	High	High	High	Croplands, Marginal lands
Biodiesel (FAME*)	Oil seeds (rape, soy)	US, Europe	Moderate	Moderate	Low	Croplands
Biodiesel (FAME)	Palm Oil	Southeast Asia	Moderate	Low-moderate	Moderate	Coastal lands
Biodiesel (BTL**)	Biomass	None	High	High	High	Croplands, marginal lands

Source: IEA

*Fatty Acid Methyl Ester

**Biomass-to-liquids

Table 3-3: Worldwide - Top Biofuel Producing Countries (in Million Liter), 2005

	Ethanol	Biodiesel
United States	16,214	284
Brazil	16,067	70
China	3,800	-
India	1,700	-
France	910	557
Russia	750	-
South Africa	390	-
Spain	376	84
Germany	350	1921
Thailand	300	-
United Kingdom	290	74
Ukraine	245	-
Canada	230	-
Poland	220	-
Indonesia	170	-
Saudi Arabia	170	-
Argentina	165	-
Italy	150	227
Australia	125	57
Japan	113	-
Czech Republic	-	136
Denmark	-	80
Poland	-	80
Sweden	-	7
Other Countries	2,139	463
World	44,874	4125

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Renewable Energy in Brazil

- In 2007: 42% of primary energy came from renewable resources
 - 13% was Hydropower
 - 29% was Biomass
 - 13% Sugarcane
 - 13% Wood and related
 - 3% different other sources
- 50% of transportation fuel comes from Ethanol
- about 2% of arable land is used for production of bioethanol from sugarcane

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
Field numbers

Use of land in Brazil

The crop of sugar cane is growing in Brazil, but according to Conab, it is not a threat to the area where food is produced

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Bioethanol

Global Availability WORLD MAP OF SUGARCANE PRODUCTION

100 countries could supply BioFuels to 200 nations,
while currently 20 oil producers provide fossil fuels to 280 nations.

— Source: British Sugar, Prepared by UNICA —

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


Figure 3-8: Worldwide - Forecast of Ethanol Demand (in Million Cubic Meter), 2007-2010

Source: F.O. Litch, RNCOS

Figure 4-1: Brazil - Production of Vehicles (in Million), 1980, 1985, 1990, 1995, 2000 & 2005



Source: Brazilian Association of Car Manufacturers

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
Different feedstocks for Bioethanol...

The global targets for Biodiesel are...

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
Different feedstocks for Biodiesel

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Current Global Biofuel Facilities ...

	Biodiesel	Bioethanol
Production sites		
North America	165	200
South America	60	340
European Union	120	70
India	15	25
China	35	50
Rest of World	60	40
Total	455	545
Blenders	170 +	
Regulator/TestLab	650 +	
Total	~ 2000*	

*Universities and Agricultural-R&D facilities are not included; this would add another 300-500

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

Algae

Poplar

- Food vs. Fuel
- Marginal Lands
- More Efficient
- Reduce the cost to produce and distribute

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
Objectives

- To develop a method for the simultaneous, high through put and accurate analysis of trace metal impurities in biodiesel samples by using ICP-OES (EN 14538, EN 14107)
- Perform analysis at the levels specified in ASTM D 6751 and EN 14214
- Apply the developed method for the quality control of biodiesel samples obtained from different feed stocks.

Metal impurity	Effects
Phosphorus	Can damage the catalytic convertor system
Calcium & Magnesium	Contribute to injector clogging, fuel pump, piston and engine wear, engine deposits
Sodium & Potassium	Contribute to injector clogging, fuel pump, piston and engine wear, engine deposits Can cause soap formation May get collected in exhaust Adverse effects on particulate trap technologies
All metals	Will affect oxidational stability of biofuels

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
ANÁLISE DE METAIS POR ICP-OES

Parameter	
Plasma gas flow	18 L/min
Auxiliary gas flow	1.6 L/min
Nebulizer gas flow	0.5 L/min
RF power	1450 watts
Viewing height	15 mm
Plasma view	Radial
Read delay	90 seconds
Read parameters	Min 0.5 and max 2.0 s
Pump rate	1.0 mL/min
Processing	Peak area
Spray Chamber	Baffled cyclonic
Nebulizer	Low flow GemCone
Injector Alumina	1.2 mm i.d.
Sample tubing	Solvent flex (0.76 mm)
Drain tubing	Solvent flex (1.14 mm)
Quartz torch	3-slot
Sample capillary	Teflon® (1 mm)

- Two biodiesel samples based on different feed stocks were analysed
- Palm oil and Jatropha based
- Samples were diluted 1:10 with high purity kerosene
- Internal standard was mixed with the diluent
- Blank was prepared with PerkinElmer blank oil and diluted with diluent containing Internal standard
- Samples were prepared in duplicate
- QC sample was prepared from Spex Certiprep ICP oil standard
- One sample from each feedstock was spiked with QC sample prior to analysis

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U.S. and EU specifications for elemental content in biodiesel (mg/kg)

Metal	ASTM D6751	EN 14214
Phosphorus	10	10
Calcium	5 (combined)	5 (combined)
Magnesium		
Sodium	5 (combined)	5 (combined)
Potassium		

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


FIELD APPLICATION REPORT

ICP OPTICAL EMISSION

Sulphur and Phosphorus Analysis in Vegetable Oil and Beef Tallow for Biodiesel Production Using the Optima Inductively Coupled Plasma-Optical Emission Spectrometer

APPLICATION NOTE

ICP-Optical Emission Spectroscopy

Author

Praveen Sarojam, Ph.D.
Global Application Laboratory
PerkinElmer, Inc.
710 Bridgeport Avenue
Shelton, CT 06484 USA

ecoanalytix™

APPLICATION NOTE

Biofuels

Authors

Zoe A. Grosser, Ph.D.

Lee J. Davidowski, Ph.D.

Pamela Wee

PerkinElmer, Inc.
Bridgeport Avenue
Shelton, CT USA

FIELD APPLICATION REPORT

ICP-OPTICAL EMISSION SPECTROMETRY

Phosphorus, Calcium and Magnesium Analysis of Soybean Oil-Feedstock for Biodiesel Production Using the Optima Inductively Coupled Plasma-Optical Emission Spectrometer (ICP-OES)

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
➤ Weigh 100 mg +/- 0.1 mg of biodiesel into a 10 mL septa vial

Internal Standard Addition (ISTD)

- Add 100 µL of 1,4 Butantriol & Tricaprin by syringe to the vial
- Add 100 µL of MSTFA (derivitizing solution)
- Shake solution and allow to react to for 15-20 minutes

Dilute with 8 mL of Heptane and shake to mix



This solution is now ready for analyses!

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
Components

Identification | Calibration | User Values/LIMS

Component Type
☐ Peak ☒ Named group ☐ Timed group

Name:

Group Members

Glycerin
1,2,4 Butanetriol
Mono1
Total Mono
Mono2
Mono3
Mono4
Mono5
Tricaprin
Total Diglycerides
Total Triglycerides

Internal Standard:

Next Previous New Component Delete Component

[Edit Component](#) OK Cancel Apply


The name by which the component is to be identified and labeled

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
Components

Identification | Calibration | User Values/LIMS

Component Type
☐ Peak ☐ Named group ☒ Timed group

Name:

Start time: min

End time: min

☐ Include only unidentified peaks

Reference:

Internal Standard:

Next Previous New Component Delete Component

Edit Component OK Cancel Apply

The name by which the component is to be identified and labeled

Component List:
Glycerin
1,2,4 Butanetriol
Mono1
Total Mono
Mono2
Mono3
Mono4
Mono5
Tricaprin
Total Diglycerides
Total Triglycerides

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
The Mono-, Di- and Tri- glycerides have a portion of the molecule due to glycerin structure. This is the bound glycerin. Only the glycerin portion of the glycerides are to be included in the total glycerin in the biodiesel. This requires each of the calculated glycerides to be multiplied by a factor relating to the glycerin portion prior to summing all the results.

Example, per the ASTM D-6584:

Glycerin $0.015 * 1.0000 = 0.015$
Total Monoglycerides $0.016 * 0.2591 = 0.004$
Total Diglycerides $0.194 * 0.1488 = 0.029$
Total Triglycerides $0.032 * 0.1044 = 0.003$

Total Glycerin = 0.051 % mass

Make sure your data-handling package can do this without exporting to some other second party software

Components

Identification Calibration User Values/LIMS

User label: []

User Values:

1	0.00000
2	0.00000
3	0.25910
4	0.00000
5	0.00000

☐ Send results for this component to LIMS

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

Mét. EN14110 - Metanol

- Methanol is used in the reaction process with a catalyst to create the fatty acid methyl ester, glycerin is formed as a by product
- Residual methanol and glycerin are removed by washing with water 4 to 6 times or more
- Methanol can be tested by:
 - Flashpoint
 - Quantification is Highly subjective and slow
 - Infrared Spectroscopy
 - Below 1% methanol is difficult to quantitate
 - Gas Chromatography
 - EN 14110 – Headspace GC
 - ASTM methods yet to be created

Europe utilizes a HS-GC technique for determination of residual MeOH, this is a simple technique. No ASTM method exists yet. Most users follow EN 14110 or modify EN 14110.

- 5 g of biodiesel is added to a 20 mL HS vial, with 5 μ L of 2-propanol (ISTD), or ESTD procedure is also accepted
- The vial is heated and equilibrated, by the HS instrument, at 80 °C for 45 min
- Manual headspace injection is allowed but automated HS injection is recommended in the EN 14110 method

Any of several columns can be used, none are compatible in the same GC oven as the glycerin column due to the high temperature of the glycerin analysis, therefore a second GC is required.

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
- Simplest and least expensive method, but not recommended
- Utilizes same GC and column as glycerin analysis
- Biodiesel sample is diluted with heptane or injected neat
- Internal Standard can be used
- Not as accurate and slower than headspace
- Methanol is unretained and therefore co-elutes with any very minor unretained artifacts (biggest reason for not recommended)
- Probably not acceptable once ASTM method is created

Instrument Conditions:

- POC inlet: 50 °C initial, hold for 2 minutes, ramp to 380 °C
- FID: 380 °C
- Oven: 50 °C initial, hold for 2 minutes, ramp to 380 °C

Method takes approximately 9 minutes

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Biodiesel Glycerin & Methanol Analyzer

ASTM D6584, EN 14105, EN 14110


Automated headspace sampler optimizes productivity and simplifies sample preparation

Rapid Ramp-up - complete package with instrument, software, SOPs, consumables, application-specific installation and qualification procedures


Two Solutions in One: Dual oven to complete both GC analyses on a single instrument; no need to modify hardware setup or swap columns

Clarus® 500 GC
IrboMatrix™ Headspace

0.1% methanol in a biodiesel matrix

Analysis of a biodiesel sample for free and total glycerin.

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

The EcoAnalytix™ Biodiesel Glycerin & Methanol Analyzer

- All of the benefits from the Clarus 500 GC plus the very best from TurboMatrix HS
- 12-sample overlapping heating oven for high throughput
- Pressure balanced injection
- Very reliable and repeatable sampling

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Biodiesel FAME Determination by GC

EN 14103


Clarus® 500 GC

Profiling the source of feedstock, looking for unwanted by-products

Calculation of different parameters from FAME profile, like iodine number, degree of transesterification, etc.

Rapid Ramp-up of sampling and analysis with consumables package

Intuitive touch-screen interface allows quick set up and operation

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FIELD APPLICATION REPORT

GAS CHROMATOGRAPHY

Free and Total Glycerol in B100 Biodiesel by Gas Chromatography According to Methods EN 14105 and ASTM D6584

Introduction

With today's increasing concern for the environment and the depletion of fossil fuel resources comes a greater awareness for alternative fuels, especially for biofuels. One of the more common biofuels is biodiesel, which is a renewable fuel from natural oils such as soybean oil, rapeseed oil or animal fats; it is a substitute for petroleum-diesel fuel.

Biodiesel consists of fatty acid alkyl esters produced by the transesterification reaction of vegetable oils and animal fats

Authors

Tim Ruppel
Gerald Hall
Timon Huybrighs
William Goodman

PerkinElmer, Inc.
710 Bridgeport Avenue
Shelton, CT USA

APPLICATION NOTE
GAS CHROMATOGRAPHY

Free and Total Glycerin in B100 Biodiesel by Gas Chromatography

Authors:
Timothy Ruppel and Gerald Hall
PerkinElmer Life and Analytical Sciences
710 Bridgeport Avenue
Shelton, CT 06484

Fatty Acid Methyl Esters in B100 Biodiesel by Gas Chromatography (Modified EN 14103)

Introduction

The production and consumption of biofuels continues to increase as more attention is paid to the environment and the depletion of fossil-fuel resources. Biodiesel, a fuel from natural oils such as soybean oil, rapeseed oil or animal fats, is a substitute for petroleum-diesel fuel. The quality criteria for the production of biodiesel are specified in EN 14214.

Authors

Timothy Ruppel
Timon Huybrighs

PerkinElmer, Inc.
710 Bridgeport Avenue
Shelton, CT USA




APPLICATION NOTE
GAS CHROMATOGRAPHY

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Biodiesel IR FAME Analyzer


ASTM D7371, EN 14078

Wide range of accessories for many different sample types

Spectrum 100 FTIR

Rapid Ramp-up
SOPs, consumables package with starter calibration model

Typical Spectrum for B5 Blend


Identification and quality of incoming feed stocks; monitoring of conversion process to determine progress of reactions & yield

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Biodiesel blend analysis by FTIR

- **EN 14078**
 - European standard
 - Transmission (0.5 mm)
 - Linear correlation against carbonyl band

- **ASTM D 7371**
 - American standard
 - ATR sampling
 - Chemometric calibration (PLS) using two spectral regions (shaded)

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Why use ATR?

- ATR effectively allows a very short pathlength (a few μm)
- No need for sample dilution
- Sample introduction is much easier (no syringes needed)
- Cleaning is very quick and requires minimal solvent
- Diamond ATR maximises the advantages:
 - Extremely durable
 - Excellent reproducibility
 - Extremely easy cleaning

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
FAME measurement procedure

1. Place one drop of the sample on the ATR crystal
2. Press the Go button to measure the spectrum
3. Clean the crystal
4. Click the Biodiesel button in the Spectrum Express software

Samples View 1		Results Table		
		Results		
 Sample Name	Description			FAME Vol
▶ UATR-10				16.02
UATR-11				18.06
UATR-32				16.05
UATR-33				18.10
UATR-34				20.01
UATR-54				15.87
UATR-55				18.15
UATR-56				20.12

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
Reporting

- Reports can be automatically printed or saved (.rtf format)
- A running record of analysis results can be saved to a .csv file that can be opened in Excel

PerkinElmer Spectrum Express Version 1.03.00
28/07/2008 18:52

FAME Analysis by ASTM 7371-07 Report

Analyst	Administrator
Date/Time	28 July 2008 18:52
Software Product	PerkinElmer Spectrum Express Version 1.03.00
Report File	C:\pel_data\biodiesel\ASTM_7371\Reports\2008_7_28\UATR-1_FAME_ASTM_7371(1).rtf
Sample Description	U-FAME
Low Concentration (0 - 10%) Quant Method	C:\Quant\PLS1\Methods\BDUATR.MD
Low Concentration Component	6.00668
Medium Concentration (10 - 30%) Quant Method	C:\Quant\PLS1\Methods\BDUATR.MD
Medium Concentration Component	6.00668
High Concentration (30 - 100%) Quant Method	C:\Quant\PLS1\Methods\BDUATR.MD
High Concentration Component	6.00668
Pathlength	1


Sample	RMS Value	Peak to Peak Value	Total M-Distance	Residual Ratio	FAME Value
UATR-1	0.0021852	0.0082682	0.65742	1.635	-0.28318

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“Starter” calibrations for ASTM D 7371

- Calibration models developed according to ASTM D 7371 methodology, using soy-based biodiesel
- Separate PLS models for 0–10%, 10–30%, 30–100% FAME
- Start using the analyzer immediately

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Biodiesel Oxidation Stability Testing by DSC

ASTM D6186


Unique Fast Scan DSC technology for unmatched sensitivity and superior signal resolution so you don't miss transitions

DSC 8000/8500

Quickly screen for oxidation stability in less than 10 minutes

Low-mass furnaces for quick response time and high throughput; Power compensation design for true isothermal measurements

Sample	1	2	3	4
Oxidation onset temperature (degrees Celsius)	155	167	155	134

Table 1. OT summary
SOURCE: PERKINELMER

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
Diesel Oil on the Clarus 600 with High Power Heater

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
From the Restek Catalog

■ Diesel
Rtx®-1
(split injection)

30m, 0.32mm ID, 1.5 μ m Rtx-1 (cat.# 10169)
1.0 μ L split injection of a diesel fuel #2 standard.
On-column concentration 165ng
mix cat.# 31093


Oven temp.: 40°C (hold 5 min.) to 300°C @
4°C/min. (hold 5 min.).
Inj. / det. temp.: 275°C / 300°C
Carrier gas: helium
Linear velocity: 30cm/sec. set @ 40°C
Split ratio: 30:1

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

ASTM D-2887 C5 to C44 Standard Rtx®-2887

(direct injection)

10m, 0.53mm ID, 2.65µm Rtx-2887 (cat.# 10199)
1.0µL direct injection of a 0.1 to 0.01 wt % hydrocarbon
standard in carbon disulfide.

Oven temp.: 35°C to 360°C @ 15°C/min. (hold 5 min.)
Inj. & det. temp.: 360°C
Carrier gas: nitrogen
Linear velocity: 112cm/sec. (15mL/min.)


■ Sim Dist

ASTM D-2887 Reference Gas Oil #1 Rtx®-2887

(direct injection)

10m, 0.53mm ID, 2.65µm Rtx-2887 (cat.# 10199)
1.0µL direct injection of a 0.1 to 0.01 wt % hydrocarbon
standard in carbon disulfide.

Oven temp.: 35°C to 360°C @ 15°C/min. (hold 5 min.)
Inj. & det. temp.: 360°C
Carrier gas: nitrogen
Linear velocity: 112cm/sec. (15mL/min.)

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Figure 1. The switching and splitting wafers

Figure 2. Installing the wafer in the holder

Figure 3. Wafer mounted in the threaded holder


Figure 4. The holder containing the wafer mounted on the GC oven wall and external connections being made.

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
S2 **S3** **S4** **S5** **S6**

Help Exit


This utility software is made available solely as a reference tool. It is provided only on an "AS IS" basis and on the express condition that it is neither warranted nor supported by PerkinElmer, Inc. PerkinElmer assumes no responsibility or liability for the use of this utility software.

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Makes backflushing easy – get rid of unwanted components


Significantly reduces runtime

Less wear on column; therefore, they last longer ... reduces costs

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
* not specified yet, can be measured by GC

** Intended to be measured by IOP-OES

*** not specified yet, can be measured by LC

**** not specified yet, can be measured by LSC

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Bioethanol Alcohol Determination by GC


ASTM D5501, DIN 51627-3

Fastest heat-up and cool-down
conventional GC oven optimizes
productivity


Rapid Ramp-up - consumables
package with calibration
standards provided

Clarus® 600 GC

For determination of methanol,
ethanol, other alcohols ranging
from C1 to C5 (extending to C10)
and other denaturants
(hydrocarbons and aldehydes)

Overlay of Calibration standards

Typical Chromatogram of Bioethanol

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Ethanol for Fuel - ASTM D5501

- Clarus 500/600 GC with Liquid Autosampler and programmable pressure control (PPC)
- TotalChrom control and data handling software
- Capillary Injector
- Flame Ionization Detector (FID)
 - Ethanol 93-97%
 - Methanol 0.1-0.6%
 - Higher boiling components 3-6% total

Inside the oven:

Column- Elite 1, 100-150 meter, 0.5 um film thickness, 0.25 mm ID

- Any column allowed that can give separation
 - Long columns = long analysis time
 - Some people use 30m BAC1 type column for 10 min analysis
-
- Ethanol can be blended into gasoline E10 - E85
 - Pure ethanol at Indy 500

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
Ethanol fuel contaminant compounds

	Sys1				
methanol	7.28				
acetaldehyde	7.28				
acetol (1,1-diethoxy ethane)	7.28				
ethanol 0.5min wide	8.06				
acetone	8.80				
isopropyl alcohol	9.16				
ethyl formate	9.53				
tert-butyl alcohol	9.91				
n-propyl alcohol	11.12				
2-butanol	13.02				
ethyl acetate	13.51				
isobutanol (2-methyl-1-propanol)	14.19				
1-butanol	15.71				
isoamyl alcohol (3-methyl-1-butanol)	18.70				
n-amyl alcohol (1-pentanol)	19.93				
active amyl alcohol (sec-butanol)(2-methyl-1-butanol)					
System 1: 100m x 250u x 0.5u RTx1 PONA	No resolution acetaldehyde/methanol				
25cm/sec H2 @ 40degC	Poor resolution high conc ethanol/acetone/ipoh/ethyl formate				
40(12)20deg/min250(10) =35min	marginal resolution high conc ethanol/t butanol				
0.5ul inj split 40ml/min	capable of 10ppm, maybe down to 2ppm				

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
Bioethanol Fermentation Determination by LC

Fermentation Broth Analysis

Flexar HPLC
Refractive Index Detector

Excellent reporting capabilities with
TotalChrom

Fastest method for determining ethanol
produced and monitoring of fermentation
broth – now only 10 minutes

Typical HPLC chromatogram and report

Determines wanted and unwanted by-
products during fermentation process

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
SISTEMA FLEXAR

Solvent Degasser and Mobile Phase Bottles

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
Fermentation Analysis using Aminex Column

- Fast easy assay for eight Fermentation analytes

- dextrin
- maltotriose
- maltose
- d-glucose
- l-lactic acid
- glycerol
- acetic acid
- ethanol

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
All analyses were performed on a PerkinElmer Optima 7300DV atomic emission spectrometer utilizing the axial plasma view for best sensitivity. The system was operated in high resolution mode to achieve complete separation between closely spaced wavelengths of copper and phosphorus.

Sample introduction system comprised a PFA microspray concentric nebulizer equipped with a 100 uL/min self-aspirating capillary (ESI, Omaha, NE), a quartz cyclonic spray chamber, 1.2 mm i.d. ceramic injector and 3-slot quartz torch (PerkinElmer, Shelton, CT).

Table 2. System Operating Parameters

Plasma Gas Flow	20 L/min
Auxiliary Flow	0.5 L/min
Nebulizer Gas Flow	0.60 L/min
RF Power	1500 Watts
View	Axial
Resolution	High

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
Figure 1. A large volume cyclonic spray chamber was used. Due to the "total consumption" nature of the system, no drain was required.

FIELD APPLICATION REPORT

ICP-Optical Emission Spectroscopy

Author

Daniel H. Jones
PerkinElmer, Inc.
710 Bridgeport Avenue
Shelton, CT 06484

Contributor

Nathan Sactvet
Elemental Scientific, Inc.
2440 Cuming Street
Omaha, NE 68131

A Novel Analysis of Highly Volatile Solvents using the PerkinElmer Optima 7300 DV with a Total Consumption Sample Introduction System

Introduction

The analysis of volatile solvents for trace metal contamination is important for a variety of reasons. Inorganic contaminants can damage the refinery process and certain catalysts can be poisoned making the chemistry of the process inefficient. These contaminants may also be harmful to the products the solvents are used in. They may cause premature breakdown or wear of mechanical equipment. Perhaps most importantly, these contaminants may be discharged into the environment causing harmful effects to ecosystems.

Table 2. Samples were spiked by weight with 160 ppb organo-metallic multi-element standard. Samples were then analyzed against a calibration curve made in heptane.

Analyte	Butane (1:5)		Naptha		Gasoline		Hexane		Hexadecane (1:5)	
	Conc. ppb	Spike % Recovery	Conc. ppb	Spike % Recovery	Conc. ppb	Spike % Recovery	Conc. ppb	Spike % Recovery	Conc. ppb	Spike % Recovery
Al 308.215	<50	102	<10	97	<10	95	<10	103	<50	96
Ca 317.933	<5	103	1	100	1	89	4	99	30	99
Cd 228.802	<5	104	<1	98	<1	91	<1	102	<5	98
Cr 267.716	<5	105	1	100	<1	99	1	102	5	101
Cu 327.393	<5	103	2	99	8	97	1	103	10	101
Fe 238.204	<5	103	1	99	11	94	<1	101	5	99
Mg 285.213	<5	102	2	98	<1	99	2	102	10	97
Mn 257.610	<5	104	1	99	<1	98	1	102	5	100
Mo 202.031	<15	102	<3	99	<3	84	2	101	10	97
Pb 220.353	<75	102	<1	94	11	91	16	86	15	94
Ti 334.940	<5	104	<1	100	<1	91	<1	104	<5	102
V 290.880	<5	105	<1	100	<1	88	<1	103	5	101
Zn 206.200	<5	102	23	101	29	92	<1	98	<5	97

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
Sample Introduction System

Petrochemical/Organic Solvents High Efficiency Introduction System

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
From methanol to ethanol

For many years, 100 percent methanol was used in the high-performance race car engines. In 2006, the IndyCar Series introduced a blend of 90 percent methanol and 10 percent ethanol as a transition fuel before full incorporation of ethanol. The IndyCar Series' conversion to the renewable fuel was introduced by the late Paul Dana, a driver who was killed in 2006 during a practice session, and was supported by the Ethanol Promotion and Information Council (EPIC), a non-profit alliance of ethanol industry leaders.

The use of the Clarus GC for fuel analysis

PerkinElmer tests use its Clarus® chromatography (GC) controlled by a PerkinElmer TotalChrom® chromatography data system for collecting, processing and reporting data. The Clarus GC separates the fuel to identify additives that may give one car a competitive advantage over others. The testing takes approximately 5 minutes and can detect impurities down to concentrations of 0.10 percent.

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Version 1.00.09 - Copyright 2009